

# Quotation Form

**1. Product Name: W/ FDA! LCD Digital Infrared Thermometer Non-Contact Temperature Ear Forehead**

**2. Product Photo:**



**3. Price Terms** Price is given according to the number of orders, order-delivery date and FOB-CIF conditions.

**4.Payment terms:** Wire Transfer,LC,PayPal- Visa- Master Card- Amex -Discover

**3.Delivery time:** within 7 days after received the payment.

**4.Supplier's Name and Address :** 合肥晨星保健用品有限公司

Baoh Industrial Estate,Baoh Zone,Hefei,China

GLOBAL TRADE & GulfUS Capital LLC

One Liberty Plaza,165 Broadway 23rd Floor, New York, NY 10006 ,

DUNS- [355998415](https://www.duns.com/355998415), NAICS CODE: 425120 Wholesale Trade Agents and Brokers

---

## Item specifics

Condition:	New: A brand-new, unused, unopened, undamaged item in its original packaging	Temperature Range:	32.0-43°C/90-109.4°F
Non-Contact Thermometer Type:	Infrared Thermometer	Display:	Digital
Features:	FDA APPROVED	Material:	ABS plastic electronic component
MPN:	Does Not Apply	Battery:	2 x AA alkaline, the product does not include a b
Brand:	ADP Health	Distance to Spot Ratio:	1 to 15cm(within 5cm is the best)
EST. Delivery time::	7-12 days	Country/Region of Manufacture:	United States
UPC:	Does not apply		

Backlight LCD digital display, lightweight and portable for use.

Built-in pointer, precise non-contact infrared measurement.

2 temperature units for conversion: Fahrenheit or Celsius, long press ""MODE"" more than 5s.

Measurement modes selection for various needs conveniently.

Body temperature measurement: aiming towards the forehead with the distance of 1-15cm.

Different backlight color for different temperature range under BODY mode.

Object temperature measurement: measuring the surface temperature of an object, baby bath water temperature or milk temperature, etc.

Automatic selection range and display resolution 0.1°C/0.1°F.

Automatic data hold function and auto power off function.

Selection for buzzer; set alarm value(the alarm default value is 38°C/100.4°F).

Memorization of the last 32 groups of measurement data.

**Functions Description:**

F1: when the body temperature exceeds alarm temperature(default value is 38°C), will display red backlight and alarm;

F2: calibration(default 0), you can set when there is a temperature deviation;

F3: alarm switch(ON-OFF), when turn on the alarm, the screen will appear loudspeaker icon;

F4: modes selection(default BODY), the screen will display ""Body temp"" under BODY mode.

Backlight Color Description(ONLY BODY MODE):

<32°C: red backlight: appear ""LO"" symbol;

32-37.3°C: green backlight: body temperature is normal;

37.3-F1: yellow backlight: body temperature is high but not exceed the alarm temperature;

>F1: red backlight: body temperature exceeds F1 and alarms.

**Item Name: Non-contact IR Thermometer**

Material: ABS Plastic

Display Resolution: 0.1°C(0.1°F)

Body Temperature Measurement: 32.0-43°C/90-109.4°F

Ambient Temperature Measurement: 0-60°C/32-140°F

Accuracy: (+/-)0.2°C(0.6°C error is allowed due to the infrared test and environment effect)

Measurement Distance: 1 to 15cm(within 5cm is the best)

Auto-shutdown Time: About 7s

Power Supply: 2 x 1.5V AA Batteries(Not Included)

Operating Temperature: 10-40°C/50-104°F

Storage Temperature: 0-50°C/32-122°F

Relative Humidity: ≤85%

Item Size: 95mm x 42.5mm x 156mm/3.74" x 1.67" x 6.14" (Approx.)

**Package Includes:**

1 x Non-contact IR Thermometer (Batteries are no included),

1 x User Manual (English)

**PAYMENT TERM:**

**PayPal CREDIT**



PayPal Visa Master Card Amex Discover

**RETURN: 30 day returns. Buyer pays for return shipping |**



### Features And Uses



-  In 1 Second
-  Smart Chip
-  High temperature alarm
-  32 Units Memory array
-  Non-contact infrared
-  Large screen display



Community



Family



Enterprises and Institutions

**For Order: WhatsApp- PHONE: +1973-938-2039 [globaltradenj@gmail.com](mailto:globaltradenj@gmail.com)**

**Detail info: <http://gulfuscapital.com/medical-products/>**